

# Sorvall RC-6+ High Speed Centrifuge



## **APPLICATION:**

- High speed centrifugation of samples
- Processing of up to 6L of sample in 1L volumes in a single run
- Fiberlite rotors for ergonomic ease

## **POINT OF CONTACT and TRAINING:**

**On the job training (OJT) is required prior to use of this equipment.**

**Please see Dominique Joyner or Julian Forney for OJT.**

- Achieves higher speeds—up to 22,000rpm/55,200xg for faster separations
- Exceptional leading acceleration and deceleration rates deliver significant time savings and increase throughput
- Rotors featuring Bioseal technology for microbiological containment
- Powerful refrigeration system can cool all rotors to 4°C at top speeds, keeping samples from damage due to heat or freezing.
- Compatible with many Sorvall RC-6 and RC-5 series centrifuge rotors