

# Sorvall Legend Micro 21 Microcentrifuge Series



## **APPLICATION:**

- Benchtop centrifugation of samples
- Process up to 24 x 1.5/2.0mL samples

## **POINT OF CONTACT and TRAINING:**

**On the job training (OJT) is required prior to use of this equipment.**

**Please see Dominique Joyner or Julian Forney for OJT.**

- Acceleration up to 21,000 x g in 12 seconds
- Refrigeration offers cooling in just 9 minutes
- ClickSeal biocontainment rotor lid "clicks" in place, ensuring securely contained samples
- New PPSU lids have improved strength and chemical resistance